

Terazosin hydrochloride dihydrate

Catalog No: #S6701


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Terazosin hydrochloride dihydrate
Brief Description	Inhibitors
Purification	/
Target Name	Adrenergic Receptor antagonist
Calculated MW	459.92
Formulation	C ₁₉ H ₂₅ N ₅ O ₄ .HCl.2H ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: O.O.Cl.COC1=CC2=C(C=C1OC)C(N)=NC(=N2)N1CCN(CC1)C(=O)C1CCCO1

|c:4,6,12,14,t:2|Pathways: Neuroscience Receptor: α1-adrenoceptor Boling pt: Melting pt: Solubility: DMSO: 26 mg/mL (56.53 mM) Ethanol:

Note: This product is for in vitro research use only