

Tolcapone


Catalog No: #S6708

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tolcapone
Brief Description	Inhibitors
Purification	99.49%
Target Name	Transferase inhibitor
Calculated MW	273.24
Formulation	C ₁₄ H ₁₁ NO ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CC1=CC=C(C=C1)C(=O)C1=CC(=C(O)C(O)=C1)[N+](=O)[O-]
Pathways: Metabolism
Receptor: COMT
Boiling pt: Melting pt: Solubility: DMSO: 55 mg/mL (201.28 mM) Ethanol: 55 mg/mL (201.28 mM) Water:

References

1. Zrcher G, et al. Adv Neurol, 1990, 53:497-50

Note: This product is for in vitro research use only