

## Troxiptide

Catalog No: #S6710


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Troxiptide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	294.35
Formulation	C <sub>15</sub> H <sub>22</sub> N <sub>2</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: COC1=CC(=CC(OC)=C1OC)C(=O)NC1CCNC1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Ethanol: 3 mg/mL (10.19 mM); DMSO: 18 mg/mL (61.15 mM)  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Kusugami K, et al. Dig Liver Dis, 2000, 32(4), 305-311.

Note: This product is for in vitro research use only