

Tyrphostin AG 879

Catalog No: #S6712


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tyrphostin AG 879
Brief Description	Inhibitors
Purification	98.00%
Target Name	EGFR inhibitor; HER inhibitor; PDGFR inhibitor; Trk receptor inhibitor
Calculated MW	316.46
Formulation	C18H24N2OS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others
 SMILES: CC(C)(C)C1=CC(=C/C(=C(/N)\S)/C#N)C=C(C1=O)C(C)(C)C
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors; JAK/STAT signaling
 Receptor: EGFR; HER2-Neu; PDGFR; Trk
 Boling pt: Melting pt: Solubility: DMSO: 36 mg/mL (113.75 mM) Ethanol: 3 mg/mL (9.47 mM) Water:

References

1. Zhou Y, et al. Cancer Res, 2005, 65(13), 5848-5856.

Note: This product is for in vitro research use only