

Ulipristal

Catalog No: #S6713


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ulipristal
Brief Description	Inhibitors
Purification	/
Target Name	Estrogen/progestogen Receptor modulator
Calculated MW	433.58
Formulation	C30H37NO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:[H][C@@]12CC[C@](O)(C(C)=O)[C@@]1(C)C[C@H](C1=CC=C(C=C1)N(C)C)C1=C3CCC(=O)C=C3CC[C@@]21[H]

|r,c:16,18,24,30,t:14|Pathways:Endocrinology/HormonesReceptor:PRBoiling pt:Melting pt:Solubility:DMSO:83 mg/mL (191.42 mM)Ethanol:14 mg/mL (32.28 mM)Water:

Note: This product is for in vitro research use only