

Vinflunine Tartrate

Catalog No: #S6722


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vinflunine Tartrate
Brief Description	Inhibitors
Purification	/
Target Name	Microtubule Associated inhibitor
Calculated MW	967.02
Formulation	C45H54F2N4O8.xC4H6O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CC[C@@]12C=CCN3[C@@H]1[C@]4(CC3)[C@H]([C@]([C@@H]2OC(=O)C)(C(=O)OC)O)N(C5=CC(=C(C=C45))[C@]6[C[C@H]7C[C@H](CN(C7)CC8=C6NC9=CC=CC=C89)C(C)(F)F)C(=O)OC)OC.[C@@H]([C@H](C(=O)O)O)(C(=O)O)O

Pathways:Cell Cycle/CheckpointReceptor:MicrotubulesBoiling pt:Melting pt:Solubility:DMSO:100 mg/mL (103.41 mM)Ethanol:47 mg/mL warmed (48.6 mM)Water:100 mg/mL warmed (103.41 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only