

Zoledronic Acid

Catalog No: #S6739


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zoledronic Acid
Brief Description	Inhibitors
Purification	99.93%
Target Name	Rho inhibitor
Calculated MW	272.09
Formulation	C ₅ H ₁₀ N ₂ O ₇ P ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: OC(CN1C=CN=C1)(P(O)(O)=O)P(O)(O)=O
 Pathways: Cell Cycle/Checkpoint Receptor: Rho
 Boling pt: Melting pt: Solubility: 0.1M NaOH (aq): 40 mg/ml; DMSO: 0.004 mg/mL (0.01 mM)
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jagdev SP, et al. Br J Y, 2001, 84(8), 1126-1134.

Note: This product is for in vitro research use only