

## 5,5-Dimethyloxazolidine-2,4-dione

Catalog No: #S6744


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	5,5-Dimethyloxazolidine-2,4-dione
Brief Description	Inhibitors
Purification	/
Target Name	PKM inhibitor
Calculated MW	422.36
Formulation	C <sub>21</sub> H <sub>28</sub> NO <sub>3</sub> Br
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:[Br-].CC[N+](C)(CC)CCOC(=O)C(O)(C1=CC=CC=C1)C1=CC=CC=C1

[c:15,17,22,24,t:13,20]Pathways:MetabolismReceptor:PKM2Boling pt:Melting pt:Solubility:Ethanol: 84 mg/mL (198.88 mM)

DMSO: 84 mg/mL (198.88 mM)

Water: 84 mg/mL (198.88 mM)

Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only