


Description

Product Name	AS101
Brief Description	Inhibitors
Purification	/
Target Name	IL Receptor inhibitor
Calculated MW	312.05
Formulation	C ₂ H ₄ Cl ₃ O ₂ TeH ₄ N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1CO[Te](O1)(Cl)(Cl)Cl.[NH₄⁺] Pathways: Others Receptor: IL-1β converting enzyme Boiling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only