

Aurora A Inhibitor I

Catalog No: #S6767


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aurora A Inhibitor I
Brief Description	Inhibitors
Purification	/
Target Name	Aurora Kinase inhibitor
Calculated MW	588.07
Formulation	C31H31ClFN7O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CCN1CCN(CC1)C(=O)CC1=CC=C(NC2=NC(NC3=CC=C(C=C3)C(=O)NC3=CC=CC=C3Cl)=C(F)C=N2)C=C1
 |c:23,25,33,35,41,44,t:12,14,17,21,31,38|Pathways: Cell Cycle/Checkpoint; Chromatin/Epigenetic Receptor: Aurora A; Aurora B; Aurora C
 Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only