

NDUFA8 Polyclonal Antibody

Catalog No: #41207

Package Size: #41207-1 50ul #41207-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NDUFA8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	NDUFA8 Polyclonal Antibody detects endogenous levels of NDUFA8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NDUFA8.
Target Name	NDUFA8
Other Names	NDUFA8; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 8; Complex I-19kD; CI-19kD; Complex I-PGIV; CI-PGIV; NADH-ubiquinone oxidoreductase 19 kDa subunit
Accession No.	Swiss-Prot: P51970NCBI Gene ID: 4702
Uniprot	P51970
GeneID	4702;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of 293 cells using NDUFA8 Polyclonal Antibody



Note: This product is for in vitro research use only