

Benzenesulfonamide


Catalog No: #S6781

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Benzenesulfonamide
Brief Description	Inhibitors
Purification	/
Target Name	Carbonic Anhydrase inhibitor
Calculated MW	157.19
Formulation	C6H7NO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC=C(C=C1)S(=O)(=O)N Pathways: Metabolism Receptor: carbonic anhydrases Boiling pt: Melting

pt: Solubility: Ethanol: 31 mg/mL (197.21 mM)

DMSO: 31 mg/mL (197.21 mM)

Water:

Note: This product is for in vitro research use only