

Benzocaine hydrochloride


Catalog No: #S6782

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Benzocaine hydrochloride
Brief Description	Inhibitors
Purification	/
Target Name	Sodium Channel inhibitor
Calculated MW	201.65
Formulation	C ₉ H ₁₁ NO ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: Cl.CCOC(=O)C1=CC=C(N)C=C1 |c:10,t:5,7|Pathways: Membrane transporter/Ion channel Receptor: Sodium

Channel Boling pt: Melting pt: Solubility: Ethanol: 40 mg/mL (198.36 mM)

DMSO: 40 mg/mL (198.36 mM)

Water:

Note: This product is for in vitro research use only