

Bitopertin

Catalog No: #S6788


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bitopertin
Brief Description	Inhibitors
Purification	/
Target Name	GlyT inhibitor
Calculated MW	543.46
Formulation	C21H20F7N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C[C@H](OC1=CC=C(C=C1C(=O)N1CCN(CC1)C1=NC=C(C=C1)C(F)(F)F)S(C)(=O)=O)C(F)(F)F

|r,c:5,7,21,23,t:3,19|Pathways: Neuroscience Receptor: GlyT1; GlyT1 Binding pt: Melting pt: Solubility: Ethanol: 100 mg/mL (184.0 mM)

DMSO: 100 mg/mL (184.0 mM)

Water:

Note: This product is for in vitro research use only