


## Description

Product Name	BQ-123
Brief Description	Inhibitors
Purification	/
Target Name	Endothelin Receptor antagonist
Calculated MW	610.7
Formulation	C31H42N6O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:[H][C@@]12CCCN1C(=O)[C@@H](CC(O)=O)NC(=O)[C@@H](CC1=CNC3=C1C=CC=C3)NC(=O)[C@H](CC(C)C)NC(=O)[C@H](NC2=O)C(C)C |r,c:22,25,27,t:19|Pathways:GPCR/G ProteinReceptor:Endothelin A receptorBoling pt:Melting pt:Solubility:Ethanol: 15 mg/mL warmed (24.56 mM)

DMSO: 100 mg/mL warming (163.74 mM)

Water:

Note: This product is for in vitro research use only