

Chaetocin

Catalog No: #S6803


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chaetocin
Brief Description	Inhibitors
Purification	/
Target Name	Histone Methyltransferase inhibitor
Calculated MW	696.84
Formulation	C30H28N6O6S4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CN1C(=O)[C@@]23C[C@]4([C@@H](N2C(=O)[C@@]1(SS3)CO)NC5=CC=CC=C54)[C@]67C[C@]89C(=O)N([C@](C(=O)N8[C@H]6)NC1=CC=CC=C71)(SS9)CO)C Pathways: Epigenetics Receptor: dSU(VAR)3-9 ; mouse G9a; Neurospora crassa DIM5 Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only