

## Choline bitartrate


Catalog No: #S6805

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Choline bitartrate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	253.25
Formulation	C <sub>5</sub> H <sub>14</sub> NOC <sub>4</sub> H <sub>5</sub> O <sub>6</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C[N+](C)(C)CCO.C(C(C(=O)[O-])O)(C(=O)O)O Pathways: Others Receptor: Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only