

Coelenterazine

Catalog No: #S6808


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Coelenterazine |
| Brief Description | Inhibitors |
| Purification | / |
| Target Name | Others |
| Calculated MW | 423.46 |
| Formulation | C26H21N3O3 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:Others SMILES:OC1=CC=C(CC2=NC3=C(C(C4=CC=CC=C4)NC(=CN3C2=O)C2=CC=C(O)C=C2)C=C1
|c:8,13,15,19,31,34,t:1,3,6,11,26,28|Pathways:Others Receptor:Boling pt:Melting pt:Solubility:Ethanol: 1 mg/mL (2.36 mM)

DMSO: 84 mg/mL (198.36 mM)

Water:

Note: This product is for in vitro research use only