


Description

Product Name	CPI-637
Brief Description	Inhibitors
Purification	98.00%
Target Name	Epigenetic Reader Domain inhibitor
Calculated MW	386.45
Formulation	C22H22N6O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: C[C@@H]1CC(=O)NC2=CC=CC(=C2N1)C3=CC4=C(C=C3)N(N=C4C5=CN(N=C5)C)C Pathways: Chromatin/Epigenetic Receptor
:BRD4; CBP; EP300 Boiling pt: Melting pt: Solubility: Ethanol:

References

1. Taylor AM, et al. ACS Med Chem Lett. 2016, 7(5):531-6.

Note: This product is for in vitro research use only