

## NEDL2 Polyclonal Antibody

Catalog No: #41210

Package Size: #41210-1 50ul #41210-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	NEDL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	NEDL2 Polyclonal Antibody detects endogenous levels of NEDL2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NEDL2.
Target Name	NEDL2
Other Names	HECW2; KIAA1301; NEDL2; E3 ubiquitin-protein ligase HECW2; HECT; C2 and WW domain-containing protein 2; NEDD4-like E3 ubiquitin-protein ligase 2
Accession No.	Swiss-Prot: Q9P2P5NCBI Gene ID: 57520
Uniprot	Q9P2P5
GeneID	57520;
SDS-PAGE MW	140-150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of 3T3 cells using NEDL2 Polyclonal Antibody

---

Note: This product is for in vitro research use only