

Cyproheptadine HCl

Catalog No: #S6814


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cyproheptadine HCl
Brief Description	Inhibitors
Purification	/
Target Name	5-HT Receptor inhibitor
Calculated MW	323.86
Formulation	C ₂₀ H ₂₂ ClN ₃ O ₃ S ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: Cl.CN1CCC(CC1)=C1C2=C(C=CC=C2)C=CC2=C1C=CC=C2 |c:11,13,16,18,21,23,t:9|Pathways: GPCR/G Protein;
Neuroscience Receptor: HTBo1ing pt: Melting pt: Solubility: Ethanol: 20 mg/mL (61.75 mM)

DMSO: 65 mg/mL (200.7 mM)

Water:

Note: This product is for in vitro research use only