

## Deguelin

Catalog No: #S6817


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Deguelin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Akt inhibitor; PI3K inhibitor
Calculated MW	394.42
Formulation	C23H22O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CC1(C=CC2=C(O1)C=CC3=C2O[C@@H]4COC5=CC(=C(C=C5[C@@H]4C3=O)OC)OC)C Pathways: PI3K/Akt/mTOR Receptor: PI3K ; Akt

Boling pt: Melting pt: Solubility: Ethanol: 78 mg/mL (197.75 mM)

DMSO: 78 mg/mL (197.75 mM)

Water:

Note: This product is for in vitro research use only