

Desogestrel

Catalog No: #S6819


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desogestrel
Brief Description	Inhibitors
Purification	/
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	310.47
Formulation	C22H30O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: [H][C@@]12CC[C@@](O)(C#C)[C@@]1(CC)CC(=C)[C@]1([H])[C@@]3([H])CCCC=C3CC[C@@]21[H]
 |r,c:22|Pathways: Endocrinology/Hormones Receptor: ER; PR Boling pt: Melting pt: Solubility: Ethanol: 62 mg/mL (199.69 mM)

DMSO: 62 mg/mL (199.69 mM)

Water:

Note: This product is for in vitro research use only