

Dyngo-4a

Catalog No: #S6822


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dyngo-4a
Brief Description	Inhibitors
Purification	98.00%
Target Name	Dynamin inhibitor
Calculated MW	338.31
Formulation	C18H14N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC=C2C=C(C(=CC2=C1)C(=O)NN/C=C/3\C=C(C(=O)C=C3O)O)O Pathways: Cytoskeletal Signaling Receptor: Dynl (brain); Dynl (rec); DynII (rec) Boling pt: Melting pt: Solubility: Ethanol: 1 mg/mL warmed (2.95 mM) DMSO: 67 mg/mL (198.04 mM) Water:

References

1. Harper CB., Targeting membrane trafficking in infection prophylaxis: dynamin inhibitors. Trends Cell Biol. 2013 Feb;23(2):90-101.

Note: This product is for in vitro research use only