

## Epoxomicin

Catalog No: #S6830


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                   |   |
|-------------------|---|
| Product Name      | Epoxomicin                                      |
| Brief Description | Inhibitors                                      |
| Purification      | /   |
| Target Name       | Proteasome inhibitor                            |
| Calculated MW     | 554.72  |
| Formulation       | C28H50N4O7                                      |
| Storage           | 3 years -20°C powder; 2 years -80°C in solvent; |

## Images



## Product Description

## Research

Area: Others SMILES: CC[C@H](C)[C@H](NC(=O)[C@H]([C@@H](C)CC)N(C)C(C)=O)C(=O)NC(=O)[C@@H](N[C@@H](CC(C)C)C(=O)[C@@]1(C)CO1)[C@@H](C)O | Pathways: Proteases/Proteasome; Ubiquitination Receptor: 20S proteasome Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (180.27 mM)

DMSO: 100 mg/mL (180.27 mM)

Water:

Note: This product is for in vitro research use only