

Etretinate

Catalog No: #S6831


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Etretinate
Brief Description	Inhibitors
Purification	/
Target Name	Retinoid Receptor agonist
Calculated MW	354.48
Formulation	C23H30O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CCOC(=O)/C=C(\C)/C=C/C=C(\C)/C=C/C1=C(C(=C(C=C1C)OC)C)C Pathways: Metabolism Receptor: Boling pt: Melting pt: Solubility: Ethanol: 70 mg/mL (197.47 mM)

DMSO: 70 mg/mL (197.47 mM)

Water:

Note: This product is for in vitro research use only