

Fenoldopam mesylate

Catalog No: #S6835


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fenoldopam mesylate
Brief Description	Inhibitors
Purification	/
Target Name	Dopamine Receptor agonist
Calculated MW	401.86
Formulation	C16H16ClNO3CH4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CS(O)(=O)=O.Oc1cc=C(C=C1)C1CNCCC2=C(Cl)C(O)=C(O)C=C12 | c: 7,9,17,t: 5,21,24 | Pathways: Neuronal Signaling Receptor: dopamine-1 receptor Boiling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only