

NEIL3 Polyclonal Antibody

Catalog No: #41212

Package Size: #41212-1 50ul #41212-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	NEIL3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	NEIL3 Polyclonal Antibody detects endogenous levels of NEIL3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human NEIL3.
Target Name	NEIL3
Other Names	NEIL3; Endonuclease 8-like 3; DNA glycosylase FPG2; DNA glycosylase/AP lyase Neil3; Endonuclease VIII-like 3; Nei-like protein 3
Accession No.	Swiss-Prot: Q8TAT5NCBI Gene ID: 55247
Uniprot	Q8TAT5
GeneID	55247;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of COLO205 cells using NEIL3 Polyclonal Antibody

Note: This product is for in vitro research use only