

## Fumagillin

Catalog No: #S6841


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Fumagillin
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	458.54
Formulation	C26H34O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:Others SMILES:CC(=CC[C@@H]1[C@@](O1)(C)[C@H]2[C@@H]([C@@H](CC[C@]23CO3)OC(=O)/C=C/C=C/C=C/C(=O)O)OC)CPathways:Others Receptor:MetAP2 Boling pt:Melting pt:Solubility:Ethanol: 32 mg/mL warmed (69.78 mM)  
DMSO: 92 mg/mL (200.63 mM)

Water:

Note: This product is for in vitro research use only