


## Description

Product Name	G007-LK
Brief Description	Inhibitors
Purification	98.00%
Target Name	PARP inhibitor
Calculated MW	529.96
Formulation	C25H16ClN7O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CS(=O)(=O)C1=CN=C(C=C1)C2=NN=C(N2C3=CC=CC=C3Cl)/C=C/C4=NN=C(O4)C5=CC=C(C=C5)C#N Pathways: DNA

Damage Receptor: Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only