


## Description

Product Name	GS-9620
Brief Description	Inhibitors
Purification	99.92%
Target Name	TLR inhibitor
Calculated MW	410.51
Formulation	C22H30N6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CCCCOC1=NC2=C(C(=N1)N)NC(=O)CN2CC3=CC(=CC=C3)CN4CCCC4 Pathways: Immunology/Inflammation Receptor: TLR7 Binding pt: Melting pt: Solubility: Ethanol:

## References

1. Lanford RE. Gastroenterology. 2013 Feb 13. pii: S0016-5085(13)00169-8.

Note: This product is for in vitro research use only