

GSK J1


Catalog No: #S6847

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GSK J1
Brief Description	Inhibitors
Purification	/
Target Name	Histone Demethylase inhibitor
Calculated MW	389.45
Formulation	C22H23N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: OC(=O)CCNC1=CC(=NC(=N1)C1=CC=CC=N1)N1CCC2=C(CC1)C=CC=C2

|c:8,10,15,17,23,28,30,t:6,13|Pathways: Epigenetics Receptor: JMJD3 (KDM6B); UTX (KDM6A)

Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only