


Description

Product Name	GSK481
Brief Description	Inhibitors
Purification	/
Target Name	TNF inhibitor
Calculated MW	377.39
Formulation	C21H19N3O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CN1C2=CC=CC=C2OC[C@H](NC(=O)C2=NOC(CC3=CC=CC=C3)=C2)C1=O
 |r,c:4,6,22,24,26,t:2,15,20|Pathways: Apoptosis Receptor: RIP1 Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only