

NET-5 Polyclonal Antibody

Catalog No: #41214

Package Size: #41214-1 50ul #41214-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NET-5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	NET-5 Polyclonal Antibody detects endogenous levels of NET-5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NET-5.
Target Name	NET-5
Other Names	TSPAN9; NET5; Tetraspanin-9; Tspan-9; Tetraspan NET-5
Accession No.	Swiss-Prot: O75954NCBI Gene ID: 10867
Uniprot	O75954
GeneID	10867;
SDS-PAGE MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using NET-5 Polyclonal Antibody

Note: This product is for in vitro research use only