

ICI-118551


Catalog No: #S6860

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ICI-118551
Brief Description	Inhibitors
Purification	/
Target Name	Adrenergic Receptor antagonist
Calculated MW	313.86
Formulation	C17H27NO2HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:Cl.CC(C)NC(C)C(O)COC1=C2CCCC2=C(C)C=C1 |c:10,19,t:16|Pathways:GPCR/G ProteinReceptor:β1-adrenergic receptor; β2-adrenergic receptor; β3-adrenoceptorBoling pt:Melting pt:Solubility:Ethanol: 62 mg/mL (197.54 mM)
DMSO: 62 mg/mL (197.54 mM)

Water: 21 mg/mL (66.9 mM)

Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only