

Isosorbide Mononitrate

Catalog No: #S6863


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Isosorbide Mononitrate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	191.14
Formulation	C ₆ H ₉ NO ₆
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: [H][C@]12OC[C@@H](O[N+](=O)[O-])[C@@]1([H])OC[C@@H]2O |r| Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 38 mg/mL (198.8 mM)

DMSO: 38 mg/mL (198.8 mM)

Water: 38 mg/mL (198.8 mM)

Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only