

Kynurenic acid


Catalog No: #S6869

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kynurenic acid
Brief Description	Inhibitors
Purification	/
Target Name	AChR antagonist; GluR antagonist
Calculated MW	189.17
Formulation	C10H7NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC=C2C(=C1)C(=O)C=C(N2)C(=O)O Pathways: Others Receptor: $\alpha 7$ nicotinic acetylcholine receptor ; glutamate receptors Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only