

Lesinurad

Catalog No: #S6875


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Lesinurad |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | OAT inhibitor |
| Calculated MW | 404.28 |
| Formulation | C17H14BrN3O2S |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Others SMILES: OC(=O)CSC1=NN=C(Br)N1C1=C2C=CC=C(C=C1)C1CC1 | c: 12, 14, 16, 19, 21, t: 5, 7 | Pathways: Membrane transporter/Ion channel Receptor: OAT1; OAT3 Boling pt: Melting pt: Solubility: Ethanol: 54 mg/mL (133.57 mM)

DMSO: 80 mg/mL (197.88 mM)

Water:

References

1. Tan, et al. [abstract]. Arthritis Rheum 2013;65 Suppl 10 :1723

Note: This product is for in vitro research use only