


Description

Product Name	LJI308
Brief Description	Inhibitors
Purification	/
Target Name	S6 Kinase inhibitor
Calculated MW	368.38
Formulation	C21H18F2N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: OC1=C(F)C=C(C=C1F)C1=C(C=CN=C1)C1=CC=C(C=C1)N1CCOCC1

|c:1,4,6,12,14,19,21,t:10,17|Pathways: PI3K/Akt/mTOR Receptor: RSK2; RSK1 ; RSK3 Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only