


Description

Product Name	LLY-507
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histone Methyltransferase inhibitor
Calculated MW	574.76
Formulation	C36H42N6O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC1=CN(C2=CC=CC=C12)CCN3CCN(CC3)C4=CC=CC=C4C5=CC(=CC(=C5)C#N)C(=O)NCCCN6CCCC6 Pathways: Epigenetics Receptor: SMYD2 Boling pt: Melting pt: Solubility: Ethanol: 22 mg/mL (38.27 mM)

DMSO: 100 mg/mL (173.98 mM)

Water:

Note: This product is for in vitro research use only