


Description

Product Name	LY3009120
Brief Description	Inhibitors
Purification	/
Target Name	Raf inhibitor; VEGFR inhibitor
Calculated MW	424.51
Formulation	C23H29FN6O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC1=CC(=C(C=C1C2=C(N=C3C(=C2)C=NC(=N3)NC)C)NC(=O)NCCC(C)(C)C)F Pathways: MAPK Receptor: B-Raf; C-Raf ; A-raf ; A-raf
 Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only