

MK-4827(Niraparib) tosylate

Catalog No: #S6892


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | MK-4827(Niraparib) tosylate |
| Brief Description | Inhibitors |
| Purification | 99.87% |
| Target Name | PARP inhibitor |
| Calculated MW | 492.59 |
| Formulation | C19H20N4OC7H8O3S |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:OthersSMILES:CC1=CC=C(C=C1)S(O)(=O)=O.NC(=O)C1=CC=CC2=CN(N=C12)C1=CC=C(C=C1)[C@@H]1CCCNC1
|r,c:3,5,16,27,29,t:1,14,18,21,25|Pathways:Chromatin/Epigenetic; DNA Damage/DNA RepairReceptor:PARP1; PARP2Boling pt:Melting
pt:Solubility:Ethanol:

References

1. Wang L, et al. Invest New Drugs. 2012, 30(6):2113-20.

Note: This product is for in vitro research use only