

## Motolimod (VTX-2337)

Catalog No: #S6898


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Motolimod (VTX-2337)
Brief Description	Inhibitors
Purification	98.00%
Target Name	TLR agonist
Calculated MW	458.6
Formulation	C28H34N4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: CCCN(CCC)C(=O)C1=CC2=C(C=C(C=C2)C2=CC=C(C=C2)C(=O)N2CCCC2)N=C(N)C1

|c:13,15,20,22,t:9,11,18,33|Pathways: Immunology/Inflammation Receptor: TLR8 Boling pt: Melting pt: Solubility: Ethanol: 15 mg/mL (32.7 mM)

DMSO: 55 mg/mL (119.93 mM) warming

Water:

Note: This product is for in vitro research use only