

Palmitoylethanolamide


Catalog No: #S6926

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Palmitoylethanolamide
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR activator
Calculated MW	299.49
Formulation	C18H37NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CCCCCCCCCCCCCCCC(=O)NCCO Pathways: Metabolism Receptor: PPARα Boiling pt: Melting pt: Solubility: Ethanol: 59 mg/mL (197.0 mM)

DMSO: 59 mg/mL (197.0 mM)

Water:

References

1. Lo Verme J, et al. Mol Pharmacol. 2005, 67(1):15-9.

Note: This product is for in vitro research use only