

## Pamapimod (R-1503, Ro4402257)

Catalog No: #S6927


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pamapimod (R-1503, Ro4402257)
Brief Description	Inhibitors
Purification	98.00%
Target Name	p38 MAPK inhibitor
Calculated MW	406.38
Formulation	C19H20F2N4O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: CN1C2=NC(=NC=C2C=C(C1=O)OC3=C(C=C(C=C3)F)F)NC(CCO)CCO Pathways: MAPK Receptor: p38 $\alpha$ ; p38 $\beta$  Boling

pt: Melting pt: Solubility: Ethanol: 28 mg/mL (68.9 mM)

DMSO: 81 mg/mL (199.32 mM)

Water:

## References

- Hill RJ, et al. J Pharmacol Exp Ther. 2008, 327(3):610-9.

Note: This product is for in vitro research use only