

## Pantoprazole


Catalog No: #S6928

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pantoprazole
Brief Description	Inhibitors
Purification	/
Target Name	Proton pump inhibitor
Calculated MW	383.37
Formulation	C16H15F2N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: COC1=C(C(=NC=C1)CS(=O)C2=NC3=C(N2)C=C(C=C3)OC(F)F)OC Pathways: Transmembrane Transporters  
 Receptor: Boling pt: Melting pt: Solubility: Ethanol: 76 mg/mL (198.24 mM)

DMSO: 76 mg/mL (198.24 mM)

Water:

Note: This product is for in vitro research use only