

## NFATc4 Polyclonal Antibody

Catalog No: #41221

Package Size: #41221-1 50ul #41221-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NFATc4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	NFATc4 Polyclonal Antibody detects endogenous levels of NFATc4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human NFATc4 around the non-phosphorylation site of S676.
Target Name	NFATc4
Other Names	NFATC4; NFAT3; Nuclear factor of activated T-cells; cytoplasmic 4; NF-ATc4; NFATc4; T-cell transcription factor NFAT3; NF-AT3
Accession No.	Swiss-Prot: Q14934NCBI Gene ID: 4776
Uniprot	Q14934
GeneID	4776;
SDS-PAGE MW	95kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of COLO205 cells using NFATc4 Polyclonal Antibody

---

Note: This product is for in vitro research use only