

Pantoprazole sodium

Catalog No: #S6929


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pantoprazole sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	HIF inhibitor; Others inhibitor
Calculated MW	405.35
Formulation	C ₁₆ H ₁₄ F ₂ N ₃ NaO ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: [Na+].COC1=CC=NC(CS(=O)C2=NC3=C([N-]2)C=CC(OC(F)F)=C3)=C1OC

|c:4,12,16,22,24,t:2,10|Pathways: Angiogenesis; Others; Chromatin/Epigenetic Receptor: HIF-1α; proton pump

Boling pt: Melting pt: Solubility: DMSO: 81mg/ml(199.83 mM); Water: 81mg/ml(199.83 mM); Ethanol: 81mg/ml(199.83 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Shen Y, et al. Oncol Lett. 2016, 11(1):717-722.

Note: This product is for in vitro research use only