

Pazopanib HCl (GW786034 HCl)

Catalog No: #S6930


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Pazopanib HCl (GW786034 HCl) |
| Brief Description | Inhibitors |
| Purification | / |
| Target Name | FGFR inhibitor; PDGFR inhibitor; VEGFR inhibitor |
| Calculated MW | 473.98 |
| Formulation | C21H23N7O2SHCl |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Others SMILES: Cl.CN(C1=CC2=NN(C)C(C)=C2C=C1)C1=NC(NC2=CC=C(C)C(=C2)S(N)(=O)=O)=NC=C1

|c:9,12,24,30,32,t:2,4,15,19,21|Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: FGFR; PDGFR; VEGFR1; VEGFR2; VEGFR3

Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only