

Peficitinib (ASP015K, JNJ-54781532)

Catalog No: #S6933


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Peficitinib (ASP015K, JNJ-54781532) |
| Brief Description | Inhibitors |
| Purification | / |
| Target Name | JAK inhibitor |
| Calculated MW | 326.39 |
| Formulation | C18H22N4O2 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Others SMILES: NC(=O)C1=C(N[C@H]2C3CC4CC2C[C@](O)(C4)C3)C2=C(NC=C2)N=C1 | r,c:3,23,26,t:20 | Pathways: JAK/STAT
 Receptor: JAK Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only